



### CATALOG PRODUCTS COVID -19

In this difficult moment, Velas Largas , together with a pool of companies specialized in Marketing and International Procurement , wants to make its contribution to public and private companies that need to find protection devices on the market, which today are difficult to supply in terms of quality, certifications of compliance with European and American regulations (CE and FDA) for large orders necessary to deal with the protection required in the various work sectors.

We have also decided to contribute to the orders relating to the three -ply protective masks made by the PAs, Hospitals, RSA and Educational Institutions, by activating a donation system to supplement the orders placed .

### CONTACTS


<https://www.velaslargas.com>

Rua Dr. António Martins, Número 44, 2º Esq,  
Benfica, Lisboa, Portugal  
Tel.+39 3355767181  
info@velaslargas.com



# INITIATIVE COVID-19

Products Supply Emergency



DNA Virus Sampling



Oral Sampling



Nasal Sampling





## 1. STERILE SURGICAL MASK

### Product description

**Sterile surgical mask for medical use,** disposable, latex-free, hypoallergenic, with glass fibre, non-woven and fluid resistant.

### Three-ply construction

Wide face coverage given by the extension of the 3 folds that allow you to expand the mask so as to cover the entire area from the nose to the chin. With laces or elastic bands.

### Technical features

**99% bacterial filtering efficiency (BFE) at a capacity of 5 microns;**

It minimizes exposure and contact with blood and body fluids by blocking contamination by the patient on the aspiration of infected microorganisms.



### STERILE SURGICAL MASK

### MEDICAL USE

**Surgical Mask 3 ply sterile**  
Size: 175\*95mm

Plus 22% VAT if imported into Italy

**Minimum order in pieces**

**The supply includes:**

1,000,000 pcs / 20,000 Box / 500ctn  
500,000 pcs / 10,000 Box / 250ctn  
100,000 pcs / 2,000 Box / 50ctn

**Monthly production capacity**  
20/60 million monthly

**Carton box capacity**  
**Mask Box capacity**

40 box  
50 pcs

Box and Carton

**Certifications and Test Report**

CE, FDA

Documentation attached to the order

**Insurance policy**  
Procurement Partner

YES

The supply is guaranteed by an insurance policy

**Delivery times**

10/15 days

Working days

**Production inspection**  
Conformity and certification of the goods

Alibaba  
Bureau Veritas  
SGS

The production of 3-ply surgical masks is subject to quality control inspection with release of product conformity certification before the goods are shipped.

**Price: FOB (free on board)**  
Shipping costs and delivery times are excluded.

**Payment terms**  
non-negotiable

50% order  
50% of goods produced and inspected

The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.

**! PA, Hospitals and RSA Donation Program**

YES

For each order placed, 1% of protective masks will be donated in addition to the order



# INITIATIVE COVID-19

Products Supply Emergency



## Nano-Tech FFP2/N95 Masks



CE Certificate: IP200318.ZRM US33  
Stock: 100.000 pcs  
Max Output Capacity per month: 40.000.000 pcs

A newly innovated mask with 95% filterability was just placed into market early this year. It could reach European and U.S. standards equivalent to KN95 level. The most impressive function is that it could still keep with stability and high filtration efficiency after disinfection either by boiled water or ethyl alcohol. It could be reused up to 20 times. Filter element is made of nano fabric micro membrane material. Its filtering performance, blood impermeability and other critical data all reached international mask standards for medical protection purpose. Meanwhile, it is also anti-corrosive and environmental friendly.



## INITIATIVE COVID-19

Products Supply Emergency



### Nanotechnology-enabled N95 Masks Can Halt the Spread of Coronavirus



These days all of us are worried about coronavirus and try our very best to protect ourselves in as many ways as we hear may actually work; wearing N95 masks designed based on nanofiber technology is one of the ways that is reported to be highly effective.

As the Wuhan coronavirus continues to claim lives in almost every country of the world, especially in **China**, where it first originated, the officials of this country strongly urge people to wear face masks in public to stem the outbreak of this lethal virus.

But a key question many people are these days asking is “Can wearing a regular surgical face mask protect you against coronavirus?” The answer is “No.” To make sure you do not catch the virus, you should wear a more specialized face mask, known as an N95 respirator.

Respirators are a particular type of personal protective equipment designed to reduce the respiratory exposure of the wearer to dangerous substances such as toxic chemicals or infectious particles. What all types of these personal air purifiers have in common are the filters specifically tailored for use in different environments.

**Nanofibers** are the perfect filter material for use in N95 respirators, which not only do provide a very large specific surface area, but can also be functionalized with various chemically active groups to improve the efficiency of the respirator in capturing naturally occurring nanoparticles (e.g., viruses), micron-sized particles (e.g., bacteria), and man-made particles (e.g., soot from diesel exhaust).

Several methods are available for nanofiber production, such as conjugate spinning, chemical vapor deposition, phase separation, drawing, self-assembly, melt-blowing, and electrospinning; among which electrospinning is a widely used technique for fabricating air filter media, because electrospun nanofibers have an average specific surface area of 1000 times higher than microfibers, along with an average pore size of 4 to 100 times smaller than microfiber membranes. Hence, an electrospun nanofiber membrane has a remarkably higher capturing efficiency than a microfiber membrane, making a perfect filter for N95 respirators.



# INITIATIVE COVID-19

## Products Supply Emergency



It should be noted that N95 respirators are NOT resistant to oil and do NOT provide protection against oil-containing solid or liquid aerosol particulates. But, the good news is that there is not any oil involved in the spread of coronavirus, and if worn correctly, these masks can block out at least 95 percent of the airborne particles larger than 2.5 microns, including particles or droplets carrying viruses and bacteria; thus, protecting the wearer against respiratory tract infections and diseases. Another benefit of these respirators is having a non-specific service life, which means they can be used as long as they are not damaged or breathing resistances are not detected.





## 4. PROTECTIVE MASK N95 / FFP2 / KN95 WITHOUT VALVE

### Product description:

**4-layer protective mask without valve**, hypoallergenic, non-woven fabric, with activated carbon resistant and durable to fluids.

The mask shows characteristics to filter the inhaled and exhaled air.

Wide coverage of the face over the entire area from nose to chin, with elastic bands.

**FFP2 masks are indicated for health professionals or for people exposed to low-moderate risk.** Expiry and short-lived mask if placed in dusty and infected environments.

### Technical features:

#### Bacterial filtering (BFE) equal to 92-98%

They reduce exposure and contamination on respiration of infected microorganisms.

It protects nose and mouth from contamination with particles of average diameter of 5 μ (micron).



MASK FFP2 N95 KN95 WITHOUT VALVE		MEDICAL AND CIVIL USE
FFP2/N95/KN95 without valve		Plus 22% VAT if imported into Italy
Minimum order in pieces		<b>The supply includes:</b> 500.000 pcs / 625 cartons 100.000 pcs / 125 cartons
Monthly production capacity 20/60 million monthly		
Carton box capacity Mask Box capacity	800 pcs	Polybag
Certifications and Test Report	FFP2/ CE, N95/NIOSH FDA, FFP2/KN95	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of FFP2 masks is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## 5.. PROTECTIVE MASK N95 / FFP2 / KN95 WITH VALVE

### Product description:

**4-layer protective mask with valve**, hypoallergenic, non-woven fabric, with activated carbon resistant and durable to fluids.

The mask shows characteristics to filter the inhaled and exhaled air.

Wide coverage of the face over the entire area from nose to chin, with elastic bands .

**FFP2 masks are indicated for health professionals or for people exposed to low-moderate risk.** Expiry and short-lived mask if placed in dusty and infected environments.

### Technical features:

**Bacterial filtering (BFE) equal to 92-98%**

They reduce exposure and contamination on respiration of infected microorganisms.

It protects nose and mouth from contamination with particles of average diameter of 5 μ (micron).



MASK FFP2 N95 KN95 WITH VALVE			MEDICAL AND CIVIL USE
<b>Price each mask</b> FFP2/N95/KN95 with valve			Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>			<b>The supply includes:</b> 500.000 pcs / 625 cartons 100.000 pcs / 125 cartons
<b>Monthly production capacity</b> 20/60 million monthly			
<b>Carton box capacity</b> <b>Mask Box capacity</b>	800 pcs		Polybag
<b>Certifications and Test Report</b>	FFP2 /CE, N95/NIOSH FDA, FFP2/KN95		Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES		The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days		Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS		The production of FFP2 masks is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.			
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected		The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## 6. PROTECTIVE MASK FFP3 WITHOUT VALVE

### Product description:

**FFP3 6-layer protective mask**, hypoallergenic, **without valve**, non-woven fabric, with resistant and extremely durable activated carbon and protective against fluids.

The mask shows characteristics to filter the inhaled and exhaled air.

Wide coverage of the face over the entire area from nose to chin, with elastic bands.

**FFP3 masks are indicated for doctors, health workers or people exposed to high and continuous risk or in contact with infected environments or people.**

Disposable product if used in high risk environments. Expiry and short-lived mask if placed in dusty and infected environments. Ideal for those who work with asbestos, they guarantee protection against legionella (for short periods of time) and against viruses such as avian influenza, A / H1N1 influenza, SARS, tuberculosis and Coronavirus.

### Technical features:

**Filtering at 99.9% of bacterial filtering efficiency (BFE).** Reduces exposure and contact of contamination on the breathing process of infected microorganisms



MASK FFP3 WITHOUT VALVE		MEDICAL USE
<b>Price each mask FFP3 without valve</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 500.000 pcs / 625 cartons 100.000 pcs / 125 cartons
<b>Monthly production capacity</b> 1,5 million monthly		
<b>Carton box capacity</b> <b>Mask Box capacity</b>	800 pcs	Polybag
<b>Certifications and Test Report</b>	FFP2 /CE, N95/NIOSH FDA, FFP2/KN95	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/ 15 days	Working days
<b>Production inspection</b> Co nformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of FFP3 masks is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 7. PROTECTIVE MASK FFP3 WITH VALVE

### Product description:

**FFP3 6-layer protective mask**, hypoallergenic, with valve, non-woven fabric, with resistant and extremely durable activated carbon and protective against fluids.

The mask shows characteristics to filter the inhaled and exhaled air.

Wide coverage of the face over the entire area from nose to chin, with elastic bands.

**FFP3 masks are indicated for doctors, health workers or people exposed to high and continuous risk or in contact with infected environments or people.**

Disposable product if used in high risk environments. Expiry and short-lived mask if placed in dusty and infected environments. Ideal for those who work with asbestos, they guarantee protection against legionella (for short periods of time) and against viruses such as avian influenza, A / H1N1 influenza, SARS, tuberculosis and Coronavirus.

### Technical features:

**Filtering at 99.9% of bacterial filtering efficiency (BFE).** Reduces exposure and contact of contamination on the breathing process of infected microorganisms



MASK FFP3 WITH VALVE			MEDICAL USE
<b>Price each mask FFP3 with valve</b>			Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>			<b>The supply includes:</b> 500.000 pcs / 625 cartons 100.000 pcs / 125 cartons
<b>Monthly production capacity</b> 1,5 million monthly			
<b>Carton box capacity</b> <b>Mask Box capacity</b>		800 pcs	Polybag
<b>Certifications and Test Report</b>		FFP3 / N100	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner		YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>		10/15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods		Alibaba Bureau Veritas SGS	The production of FFP 3 masks is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.			
<b>Payment terms</b> non-negotiable		50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## 8. INFRARED THERMOMETER FOR INSTANT TEMPERATURE DETECTION

### Product description:

#### Infrared thermometer for temperature measurement.

Product of excellent manufacturing quality with measures that allow easy portability. The thermometer is equipped with LCD backlight.

Equipped with temperature range : 32 °C ~ 42.9 °C (body temperature mode) has an accuracy of  $\leq 0.3$  ° C and a maximum measuring distance of 5 to 10 cm.

#### Technical features:

Audible alarm that activates when the thermometer detects a low or high temperature .

**Very quick response time, just 0.5 ''.**

Storage of all temperatures detected.



INFRARED THERMOMETER			MEDICAL AND PERSONAL USE
Price each thermometer			Plus 22% VAT if imported into Italy
Minimum order in pieces			<b>The supply includes:</b> 200.000/ 10.000 cartons 50.000/2.500 cartons
Monthly production capacity			
Carton box capacity		20 pcs	Box and Cartons
Certifications and Test Report		CE, FDA	Documentation attached to the order
Insurance policy Procurement Partner		YES	The supply is guaranteed by an insurance policy
Delivery times		10/15 days	Working days
Production inspection Conformity and certification of the goods		Alibaba Bureau Veritas SGS	The production of infrared thermometer is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.			
Payment terms non-negotiable		50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 9. TOTAL PROTECTION ANTIBACTERIAL PROTECTIVE SUIT

### Product description:

Protective suit for medical use made of superior quality «non-woven» fabric antistatic and antibacterial light and easy to wear and with breathable polyethylene film (15-25 / m<sup>2</sup>) and melted glue (2 gr / m<sup>2</sup>) with hem, cuff and hood and 5 mm elastic.

**Soft, resistant and breathable fabric, with double zip for hermetic closure.**

Cuffs, elastic waistband and shoe cover offer the maximum possible motility with a guarantee of total protection.

### Technical features:

Double-sided tape width: 1 cm.  
 Length of double-sided tape: 67.5 cm- 74.5 cm  
 Zipper length from 68 to 73 cm (depending on requests).  
 Total fabric weight: 25/**55**/70gr./m<sup>2</sup>  
 (Quoted 55 gr.)



SIZE	WIDTH	HEIGHT
S	85-95 cm	164-170 cm
M	92-100 cm	168-176 cm
L	100-108 cm	174-182 cm
XL	108-116 cm	180-188 cm

TOTAL PROTECTION ANTIBACTERIAL PROTECTIVE SUIT		MEDICAL AND PROFESSIONAL USE
<b>Antibacterial protective suit</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 20 0.000/4.000 cartons 25.000/1.000 cartons
<b>Monthly production capacity</b> 5 million monthly		
<b>Carton box capacity</b>	25 pcs	Polybag
<b>Certifications and Test Report</b>	CE, FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/ 15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of protective suit is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## MEDICAL AND PERSONAL USE

Plus 22% VAT  
if imported into Italy

**The supply includes:**  
200.000/2.100 cartons (120 ml)  
50.000/520 cartons (120 ml)

**The supply includes:**  
50.000/1.380 cartons (da 500 ml)  
50.000/3.125 cartons (1 Lt.)

Polybag

Documentation attached to the order

The supply is guaranteed by an insurance policy

Working days

The production of sanitizing gel is subject to quality control inspection with release of product conformity before the goods are shipped.

The balance of the order must take place upon release of inspections in accordance with the goods produced before being shipped.

### HAND SANITIZING GEL 120 ml, 500 ml, 1 lt.

Sanitizing Gel 120 ml

Monthly production capacity  
5 million monthly

Sanitizing Gel 500 ml e 1 lt

Monthly production capacity  
5 million monthly

Carton box capacity

96 pcs 120 ml  
36 pcs 500 ml  
16 pcs da 1 lt.

Certifications and Test Report

CE, FDA

Insurance policy  
Procurement Partner

YES

Delivery times


10/ 15 gg.

**Production inspection**  
Conformity and certification of the goods  
**Price: FOB (free on board)**  
Shipping costs and delivery times are excluded.

Alibaba  
Bureau Veritas  
SGS

Payment terms  
non-negotiable

50% order  
50% of goods produced and inspected



## 10. HAND SANITIZING GEL 120 ml, 500 ml, 1 lt.

### Product description:

Liquid disinfectant gel for hands with high antibacterial efficacy equal to 99%.

The gel quickly eliminates bacteria and inhibits their regrowth.

Quick drying, it does not require washing and is easy to transport for its size of 120 ml (recommended for personal use), 500 ml or 1 Lt. (suitable for professional use).

All Gel formats of different capacities are equipped with a dispenser.

### Technical features:

#### 75% alcoholic base

Composition:  
75% min. alcohol  
22% water; 2% tocopherol  
1% glycerine.

High antibacterial efficacy.

Volatilization time: 5 "



SIZE TABLE

6-7cm	S
7-8cm	M
9-10cm	L

## 11. NITRILE GLOVE WITHOUT POWDER

### Product description:

Nitrile butadiene gloves for examination without powder, high sensitivity and blue color.

With anatomical shape and ergonomic shape, they allow maximum grip and protection from external agents.

Package consisting of ambidextrous gloves



### Technical features:

- Elongation:  $\geq 400\%$
- Primary skin irritation test (ISO 10993 -10);
- Skin sensitization test (ISO 10993 -10);
- Permeation (EN 374 -3);
- Viral penetration test (ASTM F 1671 -97b).

Sizes available from S, M, L, XL



NITRILE GLOVES WITHOUT POWDER		MEDICAL AND PERSONAL USE
<b>Nitrile Gloves</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 2.000.000/2.000 cartons 1.000.000/1.000 cartons
<b>Monthly production capacity</b> 5 million monthly		
<b>Carton box capacity</b>	10 Boxes each 100 pcs	Packaging Box 1000 pcs each carton
<b>Certifications and Test Report</b>	CE, FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	20/30 gg.	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of gloves is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





12. VINYL GLOVES

Product description:

High quality Vinyl gloves in transparent powder-free color. With anatomical shape, they allow maximum grip and protection from external agents. Package consisting of ambidextrous glove.

Technical features:

Length 270 /280 mm  
 Thickness 0.18 mm  
 Sizes available from S, M, L, XL



SIZE TABLE

6-7cm	S
7-8cm	M
9-10cm	L

VINYL GLOVES		MEDICAL AND PERSONAL USE
Vinyl Glove		Plus 22% VAT if imported into Italy
Minimum order in pieces		<b>The supply includes:</b> 2.000.000/2.000 cartons 1.000.000/1.000 cartons
Monthly production capacity	5 million monthly	
Carton box capacity	10 Boxes each 100 pcs	Packaging Box 1000 pcs each carton
Certifications and Test Report	CE, FDA	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	20/30 gg.	Working days
Production inspection Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of gloves is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 13. SAFETY GOGGLES

### Product description:

Protective glasses in polycarbonate and with anti-fog lenses on both sides, anti-scratch, anti-saliva, unbreakable, UV400, for medical use.

### Technical features:

Lens size: 143 mm \* 62 mm  
 Curve : 0C  
 Front length : 15,5 cm,  
 Rear frame length : 17,5 cm,  
 Frame height: 7 cm.  
 Frame thickness: 6 cm  
 Band length : 38 cm.  
 Frame Material: PVC  
 Frame color: transparent  
 Lens color: transparent  
 Product weight: 66 g / pc



SAFETY GOGGLES		MEDICAL USE
<b>Safety Goggles</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 25.000 /125 cartons 200.000/500 cartons
<b>Monthly production capacity</b> 1 million monthly		
<b>Carton polybag capacity</b>	200 pcs each carton	Polybag (like photo)
<b>Certifications and Test Report</b>	CE , FDA, EN166, ANSI Z87.1, CSA, AS/NZS1337	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	20/ 30 days	Working days
<b>Production inspection</b> Co nformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of goggles is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non -negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 14. PROTECTIVE FACIAL SHIELD

### Product description:

Protective face shield in 100% PP with double anti-fog function of transparent color.

Hypoallergenic foam band for sweat absorption. **Designed for the space necessary to maintain protective glasses or safety glasses, the protective face shield allows greater air flow and comfort.**

Welded band that gives greater resistance and reliability.

Available both with elastic band and with adjustable band and Velcro, the product proves to be a precious **light and comfortable device to wear.**

**Excellent protection from potential contamination from blood pathogens, body fluids or harmful chemical splashes.**

### Technical features:

Thickness: 0.18 mm

### Available sizes:

32 \* 19 cm

32 \* 22 cm

34 \* 26 cm



PROTECTIVE FACIAL SHIELD		MEDICAL AND PROFESSIONAL USE
Protective facial shield		Plus 22% VAT if imported into Italy
Minimum order in pieces		<b>The supply includes:</b> 25.000 - 125 cartons 200.000 - 500 cartons
Monthly production capacity 5 million monthly		
Carton polybag capacity	200 pcs each carton	Polybag (like foto)
Certifications and Test Report	CE, FDA	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of facial shield is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## 15. INSULATING SURGICAL GOWN

### Product description:

**Disposable insulating gown in non-woven PP + PE waterproof with elastic knitted cuff.**

Breathable gown, light to carry and easy to wear through the closure of the laces.

**Waterproofed with fluid insulation, it is mainly aimed at nurses, patients, doctors and laboratory operators.**

### Technical features:

Non-woven material SBPP, SMS, PP + PE in white, blue, green, yellow.  
Weight: 45 gr.



INSULATING SURGICAL GOWN		MEDICAL USE
<b>Insulating surgical gown</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 25.000/250 cartons 200.000/2.000 cartons 500.000/5.000 cartons
<b>Monthly production capacity</b> 5 million monthly		
<b>Carton polybag capacity</b>	100 pcs each carton	Polybag
<b>Certifications and Test Report</b>	CE , FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of insulating gown is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.

SIZE	WIDTH	HEIGHT
S	130 cm	105 cm
M	135 cm	110 cm
L	137 cm	115 cm
XL	140 cm	120 cm
XXL	145 cm	125 cm



## 16. SURGICAL HEAD CAP

### Product description:

The disposable surgical headgear is made of non-woven fabric.

Easy to adjust and comfortable to wear. Mainly used for the surgical operating room of the hospital.

### Technical features:

**Material:** polypropylene material, SPP / SMS non-woven fabric.

**Application:** widely used in hospitals, pharmacies, medical personnel, hygiene-related areas, the medical sector, including chemical laboratories, hygienic applications.



SURGICAL HEAD CAP		MEDICAL USE
<b>Surgical head cap</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 1.000.000/1.000 cartons 500.000/500 cartons
<b>Monthly production capacity</b> 2 million monthly		
<b>Carton polybag capacity</b>	1.000 pcs each carton	Packaging Box
<b>Certifications and Test Report</b>	CE, FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of head cap is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 17. COVER SHOES

### Product description:

CPE medical or professional disposable shoe covers (in pairs), weight 4 gr. for operating room, non-slip, waterproof and odourless blue or green color depending on availability.

Easy to fit, equipped with elastic for fixing.

### Technical features:

Material: 30µ thick polyethylene  
 Color : blue or green  
 Size: 15 x 41 cm  
 Weight: 4 gm / pcs  
 Packaging: 100 pieces



COVER SHOES		MEDICAL USE
<b>Cover shoes</b>		
<b>Minimum order in pieces</b>		
<b>Monthly production capacity</b>		<b>The supply includes:</b> 500.000/500 cartons 1.000.000/1.000 cartons
<b>Carton box capacity</b>	1.000 pcs each carton 100 pcs each box	Packaging Box
<b>Certifications and Test Report</b>	CE , FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of cover shoes is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





**18. PROTECTIVE SHOE COVER**

**Product description:**

CPE medical or professional shoe covers (supplied in pairs), weight 4 gr. for the operating room, waterproof, water-repellent and odour less non-slip, blue, light blue or green, depending on availability.

Easy to fit, equipped with elastic for fixing. Medical device with CE and FDA certification

**Technical features:**

Material: 30µ thick polyethylene  
 Color : blue or green  
 Size: cm 15 x 41 H. 45  
 Weight: 4 gm./piece  
 Packaging: 100 pieces



PROTECTIVE SHOE COVER		MEDICAL USE
Protective shoe cover		Plus 22% VAT if imported into Italy
Minimum order in pieces		<b>The supply includes:</b> 500.000/500 cartons 1.000.000/1.000 cartons
Monthly production capacity 3 million monthly		
Carton box capacity	1.000 pcs carton 100 pcs box	Packaging Box
Certifications and Test Report	CE , FDA	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of cover shoes is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 19. OXYMETER WITH OLED DISPLAY FOR FINGER OXYGEN MEASUREMENT

### Product description:


The fingertip oximeter with Oled display allows the measurement of oxygen saturation, heart rate, pulse rate and relative display of the plethysmography wave for SpO<sub>2</sub>.

Rapid detection of values in about 5 seconds. The instrument is equipped with an anti-interference function, compact and light to carry. The product switches off automatically, in the absence of signal generation, after 8 seconds of inactivity.

### Online tool to medical test standards.

### Technical features:

Battery charge, (AAA) use 10-20 hours.  
 0.96 "bicolour OLED display,  
 Display 4 directions and 5 colours  
 Visual and audible alarm function, indication of pulse rate sound.  
 Low energy consumption (less than 30mA).  
 Power supply: 1.5 V alkaline batteries (AAA format) \* 2



OXYMETER FINGER OXYGEN MEASUREMENT		MEDICAL AND PERSONAL USE
<b>Oximeter</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 6.000/15 cartons 18.000/45 cartons
<b>Monthly production capacity</b> 2 million monthly		
<b>Carton box capacity</b>	200 pcs carton 1 pc each box	Box with instructions
<b>Certifications and Test Report</b>	CE , FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days	Working days
<b>Production inspection</b> Co nformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of oximeter is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.

# INITIATIVE COVID-19

Products Supply Emergency



## Non-contact automatic body temperature detection



Detects human face and perform high-precision infrared human temperature acquisition

Measurement Distance: 0.5m-3m

1/2.7" CMOS Hi3516DV300 2Mega pixel, 1920\*1080

Lens: 4.5mm Board Lens

Power Supply: DC12V3A

Size: 126(W) \* 41.5(H) \*230(L)mm

Net Weight: 2000g

7-inch IPS HD display IP34 rated dust and water resistant

MTBF> 50000 H - Support 22400 face comparison library and 100,000face recognition records

Support one Wiegand input or Wiegand output, Support Onvif Support fog through, 3D noise reduction, strong light suppression Image stabilization, multiple white balance modes suitable for variousfields

Support electronic voice broadcast (normal human body temperatureor super high alarm, face recognition verification results

Support English, French, Italian, Spanish, Portuguese, German, Iranian,Israeli, Russian, Korean

HIGH ACCURACY, NON CONTACTAL, AVOID TOUCH SKIN



# INITIATIVE COVID-19

Products Supply Emergency




## C19 Thermal Imaging Helmet



### Non-inductive temperature measurement helmet

- Rapid detection indoor and outdoor
- Automatic recording of the user and his/her body temperature in real time
- Quickly check vehicles and fevered people
- Fast face recognition, identity verification
- Automatic recording of the user and his/her body temperature in real time



# INITIATIVE COVID-19

Products Supply Emergency



## *Application Modes:*

### **Single-Person Temperature Measurement Mode**

Quick temperature measurement for a single person

### **Large-crowd Temperature Measurement Mode**

Quick temperature measurement for multiple people

### **QR Code Recognition Mode**

Scan the QR code that carries personal information, and record personal information into the database automatically in real time, realizing personnel identification and paperless registration.

### **QR Code Recognition Temperature Measurement Mode**

Scan the QR code that carries personal information, and complete the personal temperature measurement at the same time, and record personal and temperature information in the database automatically in real time, realizing personnel identification and paperless registration.

### **License Plate Recognition Mode**

Recognize vehicle license plate, and give an early warning for any risk.

### **License Plate Recognition Temperature Measurement Mode**

Recognize vehicle license plate, and give an early warning for any risk, and complete the temperature measurement for the people in the car.

### **Human Detection Mode**

Perform thermal imaging detection on specific parts of the human body for qualitative analysis of the location and size of human lesions.

### **Facility Inspection Mode**

Perform thermal imaging inspection for industrial facilities such as HVAC equipment, pipelines, and electronic equipment. It is suitable for fire control, electricity, HVAC, and other fields, and it can also be used for night vision inspection scenarios such as finding objects, viewing people, and fire prevention.





## 20. SCANNER FOR FACIAL RECOGNITION AND TEMPERATURE MEASUREMENT

### Product description:

Facial recognition scanner, for the identification of employees, collaborators and guests, and instantaneous measurement of body temperature.

Installable at every access to the workplace it allows to recognize employees and visitors separately and to automatically enable the entrance and opening of the gates.

The collected data can be viewed and managed through a database that collects up to 30,000 names. Equipped with alarm system when temperatures are exceeded without fever (pre-set 37.3 °).

Suitable for airports, bus and railway stations, schools, hospitals, offices, commercial activities.

Facial recognition is able to detect the employee even with a protective mask distinguishing between employees and visitors. The program manages up to 24 access points simultaneously and allows you to create a historical access database



SCANNER FOR FACIAL RECOGNITION AND TEMPERATURE MEASUREMENT		PROFESSIONAL USE
<b>Scanner</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 5 pcs/1 carton 10 pcs/2 cartons
<b>Monthly production capacity</b> 1 million monthly		
<b>Carton box capacity</b>	1 pc a car tone	1 box
<b>Certifications and Test Report</b>	CE , FDA	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of scanner is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.






## 20.1 SCANNER FOR FACIAL RECOGNITION AND TEMPERATURE MEASUREMENT

### Technical Features:

<b>OPERATING SYSTEM</b>	ANDROID 7.1
<b>Processor</b>	Dual-Core 1,8 GHz + Quad-Core 1,4 GHz
<b>Memory</b>	4 GB DDR, 16 GB eMMC
<b>Display</b>	8 inches, 800 * 1280
<b>Key</b>	Reset button
<b>Sensor</b>	Infrared sensor
<b>Camera</b>	Front and dual lens (RGB + IR), 2MP + 1.3MP
<b>LED</b>	White clear
<b>Card reader without contact</b>	ISO14443, Type A/B, Mifare, Compliant ISO18092
<b>Communications</b>	Ethernet / Wi-Fi / Bluetooth
<b>Measurement accuracy of body temperature</b>	± 0,5 °, optimal operating distance 0,5 m-0,7 m;
<b>Peripheral ports</b>	1 gate DC, 1 RJ45,1 micro USB,1 Wiegand, 1 RS485
<b>Audio</b>	Digital audio speaker, Microphone
<b>Power supply</b>	12V / 2A
<b>Bracket (optional)</b>	Bracket for wall mounting Bracket for turnstile mounting or upright bracket
<b>Dimensions</b>	cm. 27,2 (L) * 13,5 (W) * 3(H)
<b>Peripheral ports</b>	-10 °C ~ 45 °C

Mobile device management (optional)

MDM



# INITIATIVE COVID-19

Products Supply Emergency



## Maghara - Walk Through Body Temperature Detector

MAGHARA portal that measures the temperature of the human body and the metal detector, is equipped with a non-contact IR temperature detection system, light during transport and simple operation for all users. It solves the problem of rapid and contactless screening of abnormal temperature targets, with high screening accuracy, greatly reducing the likelihood of staff being exposed to infections and improving detection and elimination efficiency in high crowds.




- \* Non-contact detection mode, detection distance 2 ~ 5cm
- \* Real-time display of body temperature degree
- \* Visual alarm for abnormal temperature over 37 ° C
- \* Abnormal temperature accounting over 37 ° C
- \* The adoption of single infrared detection, can be a dual sensor (optional)
- \* Harmless to the human body

### Technical specifications

- 1 Display LED Digital Display
- 2 Infrared Sensing Height 150cm, can detect 1.45~1.85m people
- 3 Infrared Sensing Built-in 1 Detector(optional for 2)
- 4 Temperature Display Real Time, within 1S
- 5 Alarming and Accounting  $\geq 37^{\circ}\text{C}$
- 6 Detection Position Forehead and Wrist
- 7 Temperature Detection Precision  $\pm 0.3^{\circ}\text{C}$
- 8 Passing Rate 70 Persons/minute
- 9 Alarming Zero Out Supported
- 10 Modular Design 10 minutes fast installation
- 11 Vertical Dimension 2200(H) x 870(W) x 350(D)mm
- 12 Vertical Channel Size 2000(H) x 710(W) mm
- 13 Dual Working Mode(optional) Metal Detection and Temperature Detection
- 14 Panel Material Fireproofing PVC
- 15 Certificate CE Approved, ISO14001 · ISO9001 · OHSAS18001
- 16 Carton Box Dimension Size 2300\*330\*460mm
- 17 N.W/G.W 30G/35KG
- 18 Working Environment
  - Temperature -20~45 °C
  - Relative Humidity  $\leq 90\%$ (RH)
  - Working Voltage AC100-240, 50~60 Hz
  - Location Indoor
- 19 Warranty 12 months
- 20 Application Train Station, School, Hospital, Office Building, etc.


# INITIATIVE COVID-19

## Products Supply Emergency



## Benzalkonium Chloride Aqueous Solution

This is a quaternary ammonium based biocide with an optimised alkyl chain length distribution, which can be used for the preparation of disinfectants, cleansers and timber preservation formulations, and for the production of fungicidal/algicidal washes and other similar applications.



### Chemical and Physical Characteristics

Composition: Benzalkonium chloride aqueous solution with the chain length distribution:

C12 (63-75%), C 14 (24-30%), C 16 (0-5%), MW 352.5 g/mol

Appearance: Clear, colourless to pale yellow liquid

Refractive Index (20°C): 1.4170 - 1.4300

Density (20°C): 0.975 - 0.995 g/ml

pH value (10% in water): 6.0 - 9.0

Cationic Activity: 49.0 - 51.0%

Amine value: max. 1.0 mgKOH/g

Colour (APHA): max. 100

Solubility: Fully soluble in water and low molecular weight alcohols

Stability in application\*: Stable over whole pH range and up to 120°C

Note: These characteristics do not constitute a sales specification.

### Biocidal Properties

It possesses a wide spectrum of antimicrobial efficacy including activity against the following organisms:

#### Bacteria

*Bacillus subtilis*

*Enterobacter aerogenes*

*Escherichia coli*

*Klebsiella pneumonia*

*Proteus vulgaris*

*Pseudomonas aeruginosa*

*Staphylococcus aureus*

#### Algae

*Chlorella vulgaris*

*Nostoc commune*

*Scenedesmus vacuolatus*

### Moulds/Yeasts

*Alternaria alternata*

*Aspergillus niger*

*Aureobasidium pullulans*

*Candida albicans*

*Chaetomium globosum*

*Cladosporium cladosporoides*

*Coniophora puteana*

*Coriolus versicolor*

*Epidermophyton \_occosum*

*Gloeophyllum trabeum*

*Microsporium canis*

*Penicillium funiculosum*

*Penicillium glaucum*

*Poria placenta*

*Saccharomyces cerevisiae*

*Trametes versicolor*

*Trichoderma viride*

*Trichophyton*

*mentagrophytes*

### Applications / Use Levels

It is suitable for the formulation of disinfectants, cleansers and sterilants and as a fungicidal/algicidal wash base for the maintenance or pre-decoration treatment of masonry, brick, timber, plaster, paintwork, etc. It can be used as an algicide/slimicide in swimming pools, industrial water reservoirs and cooling towers and for the preparation of wood preserving formulations intended for primary treatment and secondary preservation (EU/AWPA classes 2, 3 and 4, by dip or vacuum application as well as by spray application in closed tunnels).

The use concentration will depend on the particular application and the required degree of antimicrobial activity. More precise information regarding the concentration required for specific applications, such as timber preservation, can be advised or determined by the local Microbiological Technical Centre.

Prior to use it should not be diluted or mixed with other chemicals, in order to avoid any negative influences on the active ingredient(s).







## 22. VC70 PULMONARY VENTILATOR FOR COVID-19 INTENSIVE THERAPY PATIENTS

### Product description:

This lung ventilator belongs to one of the most modern ICU ventilation machines. Equipped with a metal-based heated exhalation valve, consisting of a built-in flow sensor. It is driven by a turbine and electronically controlled. The machine is suitable for infants and adults and offers two different uses for invasive and non-invasive ventilation.

Functions include Apnea ventilation, intelligent suction, manual breathing, Inspiration and Expiration block, ETCO2 measurement, nebulization, waveform freezing.

It is an intelligent machine that is capable of automatically adjusting itself, offering the patient the type of ventilation needed, thus increasing breathing comfort.

Complete monitoring of lung mechanics.

The tidal volume has a minimum of 20 ml. in VCV mode. It has 72 hours of trend with 30 seconds of resolution. Compensation of the tube, of the losses and of the conformity of the automatic circuit. More importantly, it has updated ventilation system software, with USB ports available.



PULMONARY VENTILATOR FOR INTENSIVE THERAPY		MEDICAL USE
<b>Pulmonary ventilator for intensive therapy</b>		Plus 22% VAT if imported into Italy
<b>Minimum order in pieces</b>		<b>The supply includes:</b> 15 pulmonary ventilator 1 pc for pallet
<b>Monthly production capacity</b> 25 0,000 pieces monthly		
<b>Pallet</b>	1 pallet pc	Pallet reinforced
<b>Certifications and Test Report</b>	CE FDA/MDSAP ISO	Documentation attached to the order
<b>Insurance policy</b> Procurement Partner	YES	The supply is guaranteed by an insurance policy
<b>Delivery times</b>	10/15 days	Working days
<b>Production inspection</b> Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of ventilator is subject to quality control inspection with release of product conformity certification before the goods are shipped.
<b>Price: FOB (free on board)</b> Shipping costs and delivery times are excluded.		
<b>Payment terms</b> non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 22.1 VC70 PULMONARY VENTILATOR FOR COVID-19 INTENSIVE THERAPY PATIENTS

### Technical features:

- Apnea ventilation
- Pulmonary Mechanics
- Trigger pressure and flow
- Automatic compensation of the tube
- Smart suction
- Manual breathing
- Exhalation and inspiration
- Waveform lock
- Nebulizer
- Non-invasive ventilation (NIV)
- EtCO<sub>2</sub> measurement

### Invasive and non-invasive ventilation modes:

- IPPV, SIMV, SIGH, SPONT / CPAP + PSV
- VCV, PRVC, SIMV, SIMV / VCV + PSV, PCV, APNEA, NIV / CPAP, BIPAP, manual & SGH.
- With BIVENT / APRV + PSV, NIV-T, NIV-S / T and spray treatment.

### Parameters

- Ventilation Frequency (f) VCV:
  - 1-80 minutes for adults
  - 2-80 minutes paediatric
- NIV frequency (f):
  - 4-20 minutes for adults
  - 4-40 minutes paediatric

### Inspiration times:

- 0.2-9 seconds for adults
- 0.2-5 seconds paediatric

### Current volume:

- 50-2000 ml for adults
- 20-300 ml paediatric

### T / Pause:

- 0-4 seconds for adults
- 0-2.5 seconds paediatrics
- Ventilation frequency:
  - 1-40 / minute

### Inspiratory pressure: 5-70 cmH<sub>2</sub>O

**Psupp:** 0-70 cmH<sub>2</sub>O

**Psupp NIV:** 0 to (50 PEEP) cmH<sub>2</sub>O

**Peep:** 0-35 cmH<sub>2</sub>O

**Peep NIV:** 2-20 cmH<sub>2</sub>O

**C-Pap in NIV:** 2-20 cmH<sub>2</sub>O

**Rise time (slope):** 0-2 seconds

**O<sub>2</sub> concentration (FiO<sub>2</sub>):** 21% -100%

### Sensitivity:

from 0.5 to 20L / minutes

### Flow tracking:

-20 to 0 cm H<sub>2</sub>O Trigger expiratory pressure and sensitivity (% of peak flow): 5% -80%

**Tube:** ET endotracheal tube. Tracheostomy tube



### Degree of compensation 0-100%

### Airway pressure measurement and values:

- Plateau pressure (P<sub>plat</sub>)
- Positive end expiratory pressure (PEEP)
- Peak inspiratory pressure (P<sub>peak</sub>)
- Average airway pressure (P<sub>mean</sub>)
- Min. Airway pressure (P<sub>min</sub>)

### Flow measurement:

- Expiratory minute volume (MVe)
- Spontaneous expiratory minute volume (MVespont)
- Inspiratory tidal volume (VTi)
- Expiratory tidal volume (VTe)

### Frequency measurement:

- Total respiratory rate (ftot)
- Spontaneous respiratory rate (fspont)

### O<sub>2</sub> measurement:

- Inspiratory CO<sub>2</sub> measurement O<sub>2</sub> (FiO<sub>2</sub>)
- Final expiratory CO<sub>2</sub> concentration (etCO<sub>2</sub>)



## 22.2 VC70 PULMONARY VENTILATOR FOR COVID-19 INTENSIVE THERAPY PATIENTS

### Curve display:

Airway pressure (t) (Paw) from -20 to 100 cm H<sub>2</sub>O flow (t) from -180 to 180 L / minute  
Volume (t) from 0 to 3000 ml.  
CO<sub>2</sub> (t) from 0 to 80 cm H<sub>2</sub>O  
Pressure-volume circuit  
Pressure-flow circuit  
Flow-volume ring

### Alarms

Expiratory minute volume (MVe)  
Airway pressure (Paw)  
Expiratory tidal volume (Vte) PEEP Insp.  
O<sub>2</sub> concentration (FiO<sub>2</sub>) (automatic)

### End - expiratory concentration of CO<sub>2</sub> (etCO<sub>2</sub>) High / Low Respiratory frequency (f) High

Apnea alarm time from 10 to 60 seconds, Off Inspiration duration High  
Oxygen sensor failure Yes  
O<sub>2</sub> low supply pressure  
CA error Yes

### Low battery

Limited battery capacity Occlusion  
CO<sub>2</sub> sensor  
Leak out of range Fan block

### Performance data control principle:

Electronically guided and electronically controlled spraying: 30min  
Enrichment of suction oxygen:  
before sputum aspiration 3 minutes  
Post aspiration of sputum 2 minutes

### Basic flow. Maximum capacity

#### Loss compensation flow: operating data

5-25 L / minutes and 180 L / minutes 60L / minutes

### Main power connection.

100 V-240 V, 50 / 60Hz Energy consumption.  
110-120Vac. 2A fan only, 200VAC J  
10A fan plus 220-240 VAC auxiliary outputs.  
1A, only 200VACJ fan, 8A fan plus auxiliary outputs.

### Gas supply:

#### O<sub>2</sub> pressure gauge 280-600KP

Languages: Chinese, English, Spanish, Polish, Turkish, Russian. Other languages can be customized.

### Physical specifications:

Dimensions L \* P \* A: Model VG70  
Trolley two sizes:  
350mm \* 460mm \* 415mm  
547mm \* 675mm \* 950mm

### Weight model VG70:

Trolley 15kg  
25kg fan

### Diagonal screen size:

12 " color TFT touch screen  
Input / output ports with 2 USB ports  
Real-time VGA output  
1 RJ 45 Ethernet connector  
Nurse call  
CO<sub>2</sub> holder







**22. MULTI-FUNCTIONAL 510S PULMONARY VENTILATOR PORTABLE FOR AMBULANCE**

**Product description:**

Pulmonary ventilator, compact, easy to carry, designed for first aid in an ambulance or hospital or for transit beds.

Methods of use and ventilation parameters can be set to facilitate the waveform of the pressure of the state of the respiratory airways.

Multiple and intuitive parameter monitoring.

**How to use:**

Meets the respiratory needs of adults, paediatric patients and infants.

**It helps optimize the quality of care and patient comfort, improving respiratory results.**

It combines advanced technology and different pneumatically controlled electronically controlled ventilation functions.



PULMONARY VENTILATOR PORTABLE FOR AMBULANCE		MEDICAL USE
Pulmonary ventilator portable		Plus 22% VAT if imported into Italy
Minimum order in pieces		<b>The supply includes:</b> 25/25 cartons re-enforced 50/50 cartons re-enforced
Monthly production capacity 10 pieces monthly		
Carton capacity	1 p.c. each carton	Carton re-enforced
Certifications and Test Report	CE FDA/MDSAP ISO	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of ventilator is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.



## 23. 1 MULTI-FUNCTIONAL 510S PULMONARY VENTILATOR PORTABLE FOR AMBULANCE

### Technical features:

With integrated battery and HD LCD display  
Pneumatic power control  
Respiratory Modes:  
A / C, A / C + Sigh, SIMV, SPONT, Manual  
Current volume: 0-1500 ml  
FiO2: 48% -100%  
P-Trig: -2-OKPa  
Mute alarm: s120s  
P-Limit -H: 2 - 6KPa  
P-Limit-L: 0-2KPa  
Peep: 0.5 - 2KPa (optional)  
Power supply: AC100 - 240 V + 10%, 50/60 HZ + 1%  
Pressure control and volume control  
Invasive and non-invasive ventilation  
SPONT / PSV, CPAP mode  
Manual intervention functions  
Manual breathing  
Inspiratory intake  
FiO2 from 40% to 100% continuously adjustable  
Display: pressure time, flow time.





## 24. PORTABLE PULMONARY VENTILATOR REHABILITATION PATIENTS

### Product description:

New BiPA P equipment with backup frequency for **non-invasive ventilation treatment**.

With 3.5-inch color LCD screen, an SD data memory card and alarms .

**Real-time visualization of pressure waves and respiratory flow.**

3-level technology that reduces pressure at the beginning of exhalation and returns to the prescribed pressure at the end of exhalation, **making breathing more comfortable and natural.**

With optional humidifier that can be connected to the central unit it forms a single integrated and complete unit to support non-invasive breathing.

Optional SPO2 module that can be connected to the BPAP unit, where it is possible to view the parameters simultaneously on the screen of the unit and in medical reports (memory card and i-Code) .



PORTABLE PULMONARY VENTILATOR REHABILITATION PATIENTS		MEDICAL USE
PORTABLE PULMONARY VENTILATOR REHABILITATION PATIENTS		Plus 22% VAT if imported into Italy
Minimum order in pieces		<b>The supply includes:</b> 15 pcs/15 cartons 30 pcs/30 cartons
Monthly production capacity 500,000 pieces monthly		
Carton capacity	1 p.c. each carton	Carton re-enforced
Certifications and Test Report	CE FDA/MDSAP ISO	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Co nformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of ventilator is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## 24.1 PORTABLE PULMONARY VENTILATOR REHABILITATION PATIENTS

### Technical features:

**Dimensions:**

170x196x118 mm.

**Dimensions with thermo humidifier:**

290x196x134 mm.

**Equipment weight:** <1.5 Kg

**Weight with humidifier:** <2,5 Kg

**Sound level:**

30 dB operating at 10 Cm H<sub>2</sub>O.

**√VNI mode:**

S / T - T - S - CPAP

**Pressure range:**

4 - 25 cm H<sub>2</sub>O

**Adjustable initial pressure:**

RESLEX (expiratory relief) at 3 levels

**Active humidifier:**

5 levels

**Minimum inspiration time:**

0.1 S - TI Max

**Maximum inspiration time:**

0.3 S - 3.0 S



# INITIATIVE COVID-19

Products Supply Emergency



## CPAP, Bi-level Ventilator



Simple and comprehensive 2-level pressure devices for the Non Invasive Ventilation (NIV) of non-dependent patients (SAS and Chronic Respiratory Failure)

- Lightweight design
- Integrated heated humidifier
- Integrated power supply
- Easy to carry
- Easy patient / machine synchronization and adaptation

It is designed to help physicians achieve effective patient adaptation and maximum therapy compliance.

The FlowSens algorithm automatically adjusts the trigger thresholds and allows the machine to instantly respond to changes in the patient's ventilatory demand. It also detects and compensates leaks avoiding auto-triggering of the device.

The "Water Control and Humid Control" features control the humidifier output and minimize condensation in the patient circuit.

### AUTO-ON

The 'Auto-on' function allows the patient to start the treatment with a single breath in the mask. The device can be started at any time by pushing the on/off button.

### "GRAPH' IT" APPLICATION

In a few clicks and at a quick glance, users can simulate the machine's settings.

### INTELLIGENT RAMP (I-RAMP)

The intelligent ramp function (I-Ramp) makes it possible to start the treatment at an adjustable, more comfortable pressure (the 'comfort pressure'), making falling asleep easier and more comfortable. It also optimizes the treatment efficacy by adapting the time in ramp to the detection of respiratory events that indicate that the patient is asleep. Once the device detects respiratory events (indicating that the patient is asleep), it will gradually and comfortably increase the pressure to reach the prescribed pressure or switch on the auto-algorithm depending on the selected mode.

# INITIATIVE COVID-19

## Products Supply Emergency



### SPECIFIC AUTO-CALIBRATION

Choose any type of mask or tubing! Versatility and flexibility were never easier thanks to our advance resistance control. This technology will ensure that the prescribed pressure is delivered at the patients interface at all times regardless of the configuration of the pneumatic pathway (mask type, diameter and length of the circuit and the presence or absence of a humidifier).

### CONNECTIVITY

Duo and Duo ST are compatible with the OxyLink module (which synchronizes oximetry data) and can be linked to the Connect telemonitoring platform.

### HUMID CONTROL+ TECHNOLOGY

The patented Humid Control + technology minimizes condensation in the patient's circuit and maintains a stable level of humidity throughout the night.

### BATTERY OPERATION

A cigarette-lighter cable is available as an optional accessory.

- MODES: CPAP, Bi-level.
- SIZE: 20 x 23 x 12 cm.
- WEIGHT: 1,5 Kg. with power supply and humidifier.
- NOISE: 27 Decibels.
- POWER INPUT: AC: 100-240 V. DC: 13 V.
- CHARACTERISTICS: auto altitude adjustment, automatic leak compensation, pressure range 4-20 cm H<sub>2</sub>O, IPAP range 3-25 cm H<sub>2</sub>O, EPAP range 3-20 cm H<sub>2</sub>O, ramp 0-45 min., back up rate 4-25 rpm., I/E ratio.
- FUNCTIONALITIES: auto calibration, FlowSens
- DATA: respiratory rate, AHI, tidal volume estimation, patient's flow, 1 year usage data, 1 year detailed data graphics.
- DATA STORAGE: memory card, micro USB cable.
- TUBE: 22 mm., 15 mm.
- HUMIDIFICATION: humid control Plus.





**INITIATIVE COVID-19**  
Products Supply Emergency



**OC70 - Humidified high flow nasal cannula oxygen ventilator**



**1.1 Performance parameter**

Parameter	Scope	Accuracy
Oxygen concentration	21%~100%	±5%
Flow	10L/min~70L/min (high flow mode) 2L/min~25L/min (low flow mode)	±5L/min (25L/min~70L/min) ±2L/min (2L/min~25L/min)
Temperature	31°C~37°C	±2°C
Humidity	2L/min~70L/min	>10mg/L
	37°C, ≤40L/min	>33mg/L
PE drop	70L/min	<2cmH <sub>2</sub> O
PI drop	70L/min	<2cmH <sub>2</sub> O
Humidifier compliance	Humidifier compliance is not affected by the water chamber	<10mL/kpa
Leakage	70L/min	<10mL/min

**1 2.Environmental parameter**

Parameter	Operating condition	Storage condition
Temperature	18°C~30°C	-10°C~60°C
Relative humidity	15%~73% (non-condensate)	15%~95% (non-condensate)
Atmospheric pressure	83~102kPa (5,600 feet above the sea level)	—————
Pressure range of high-pressure oxygen	280kPa~600kPa	—————

# INITIATIVE COVID-19

Products Supply Emergency



## 1.3 Physical & electrical parameter

Parameter	Value
Outline dimensions	340mm L x 228mm W x 162mm H
Weight	3.3kg
AC power supply	100-240VAC~50Hz/60Hz
DC power supply	24V
AC current	2.5A Max
DC current	4.2A

## 1.4 Electrical safety classification

Parameter	Value
Safety class against electric shock	Class-II
Electric shock protection degree	BF application component
Waterproof level	IP21
Mode of operation	Continuous operation
Safety class when used in combustible anesthetics mixed with air or oxygen or helium oxide	It can't be used in combustible anesthetics mixed with air, oxygen or helium oxide



**CE CFDA**



## 25. COVID-19 TEST KIT

### Product description:

**COVID-19 IgG / IgM rapid test kit (whole blood / serum / plasma) is a lateral flow immunochromatographic test.**

The diagnostics of the Covid-19 Test Kit is of the **highest quality and is among the best and the production factory is among the world's leading medical analysis centers specializing in the production of specific tests.**

The biological sample is added to the SARS-COV-2 antigen coated particles in the test cassette. The mixture then migrates upwards on the membrane by capillary action where, in the presence of IgM or IgG SARS-COV-2 complexes, it reacts respectively with anti-human IgM antibodies or anti-human IgG antibodies.

### Technical features:

- Type of sample Serum : Blood
- Serum sample volume : 20 ul
- Result Time: 5 min.
- Accuracy: 0.98
- Sensitivity: 5.0 ng/ml

Certificate  
Ministry of Health



**The test is authorized by the Ministry of Health in Italy and was carried out with advanced equipment and advanced technologies. The factory that produces it is a first-rate diagnostic centers that offers a full range of top quality health control services.**



COVID-19 TEST KIT		MEDICAL USE
COVID-19 TEST KIT		Plus VAT
Authorized Ministry of Health Italy		<b>The supply includes:</b> 50.000/25 cartons 100.000/50 cartons
Minimum order in pieces		
Monthly production capacity 3 million pieces monthly		
Carton capacity	2.000 box each carton	Box like photo
Certifications and Test Report	CE	Documentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Conformity and certification of the goods	Alibaba Bureau Veritas SGS	The production of test kit is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.





## 25.1. COVID-19 TEST KIT

### Product description:

COVID-19 IgG / IgM rapid test kit (whole blood / serum / plasma) is a lateral flow immunochromatographic test.

The diagnostics of the Covid-19 Test Kit is of the highest quality and is among the best and the production factory is among the world's leading medical analysis centers specializing in the production of specific tests.

The biological sample is added to the SARS-COV-2 antigen coated particles in the test cassette. The mixture then migrates upwards on the membrane by capillary action where, in the presence of IgM or IgG SARS-COV-2 complexes, it reacts respectively with anti-human IgM antibodies or anti-human IgG antibodies.

### Technical features:

Type of sample Serum : Blood  
 Serum sample volume : 20 ul  
 Result Time: 5 min.  
 Accuracy: 0. 98  
 Sensitivity: 5.0 ng / ml



COVID-19 TEST KIT		MEDICAL USE
COVID-19 TEST KIT		Plus VAT
Minimum order in pieces		<b>The supply includes:</b> 50.000/25 cartons 100.000/50 cartons
Monthly production capacity 2 million pieces monthly		
Carton capacity	2.0 00 box each carton	Box like photo
Certifications and Test Report	<b>CE, FDA</b>	Do cumentation attached to the order
Insurance policy Procurement Partner	YES	The supply is guaranteed by an insurance policy
Delivery times	10/15 days	Working days
Production inspection Con formity and certification of the goods	Alibaba Bureau Veritas SGS	The production of test kit is subject to quality control inspection with release of product conformity certification before the goods are shipped.
Price: FOB (free on board) Shipping costs and delivery times are excluded.		
Payment terms non-negotiable	50% order 50% of goods produced and inspected	The balance of the order must take place upon release of inspections and certifications in accordance with the goods produced and before being shipped.

# INITIATIVE COVID-19

## Products Supply Emergency



### Viral Sample Collection Kit

Certificati: Class I, # 20200086, and CE certificated



#### COMPONENTI DEL KIT

Sample Collection Kit	# 43903, 10 preps
Sample Collection Tubes	10 individual packages with biohazard bags Each including: 1 collection tube ( red-cap) with 2mL preservative 1 cotton pad
Oral Swabs	10 pieces ( sterile and individual packed )
User manual	1

#### DESCRIPTION

The Sample Collection Kit is designed for collection, transportation and short-time storage of clinical specimen total DNA/RNA of saliva, blood, serum, stool and body fluid, etc. The preservative will lyse the organism and inactivate the cells, bacteria and virus. The stored samples can be used in subsequent clinical experiments such as DNA/RNA extraction, PCR testing or sequencing.

#### STORAGE

Before use, the product can be stored at 15-30 ° C for 24 months. After use, the sample can be stored at room temperature for 1 week, or for a longer time storage, please place it at -20 ° C and below.



# INITIATIVE COVID-19

## Products Supply Emergency



### PROTOCOLS

Before starting, Read the user manual, make sure all the operations are followed as indicated.

Firstly, make sure the collection tube is sealed well and cotton pad is dry. If not, do not use and discard.

Secondary, labelling relevant sample information on the tube. Then, following the below protocols:

1. Take out of the swab from the package, and swab deep into the oral or nasal cavity or other part of the body for sampling. The specific sampling methods as follows:

a) Nasal swab: Gently insert the swab head into the nasal condyle of the nasal passage, and then slowly turn it out after staying for a while. Wipe the other nostril with another swab if necessary.

b) Pharyngeal swab: Wipe the bilateral pharyngeal tonsils and posterior pharyngeal wall with the swab.

2. Open the collection tube and immerse the swab head in the preservative, and then break the plastic stick at lower place to make the swab head sealed in the tube, discard the tail.

3. Screw up the tube cap, and shake it up-and-down several times.

4. Place the sampling tube in a biohazard bag and seal it for transportation.






# INITIATIVE COVID-19

## Products Supply Emergency



### Coronavirus Protection Cream

A SAFE AND NATURAL INGREDIENT THAT CAN PROTECT AGAINST THE CORONAVIRUS WHILE REJUVENATING YOUR SKIN



### GF2 Skin Rejuvenation

The active ingredient in GF2 is called Hypochlorous acid (HOCl) and is a naturally produced molecule also known as our immune molecule. In fact, this is the basis for our body's immune system, and it is formed in the white blood cells. It has many health benefits and contributes to the rejuvenation of the skin but also has very powerful disinfection properties. HOCl has been shown to be 80 – 100 times more effective against viruses than bleach chlorine<sup>1</sup> and it is completely natural and safe.



### The power of Hypochlorous acid

By far the safest and most effective products are those that contain Hypochlorous acid (HOCl).

The acid targets the proteins in the virus, which weakens its power to enter healthy cells. Because Hypochlorous acid is an electrically uncharged molecule it can get inside the virus and quickly kill it without damaging healthy cells. Hypochlorous acid has also been proven to reverse inflammation and “switch off” the genes responsible for inflammation, ageing and disease.



# INITIATIVE COVID-19

## Products Supply Emergency



### Is HOCl safe?

HOCl holds the highest US FDA safety classification possible. Generally Recognized As Safe (GRAS) classification is only given to a product if it is safe to inhale, ingest and apply to the skin. It is also safe to use during pregnancy and in babies.

In fact, HOCl is considered totally safe for human use. One of the primary reasons why HOCl is not used to completely replace all toxic salts and other forms of chlorine is because it is a gas and not a liquid and is very difficult to make outside of the human body.

### Not all HOCl is the same

Not all methods of manufacturing Hypochlorous acid result in medical-grade solutions. It is important to use a product that contains HOCl that has not been manufactured using salt (sodium chloride/Na). Solutions with salt in them also contain other toxic chemicals which can harm the healthy skin layer (microbiome) and cause reactions to the skin.

### Why choose this product

It owns the proprietary information to manufacture pure HOCl immune solution (without the other toxic forms of chlorine or salt present in the solution) outside of the human body. It makes HOCl that is formulated to the correct pH to achieve maximum killing power without any negative impact on the skin, eyes or mouth. It is presented in a purified water solution which makes it easy to use.



### GF2 Skin Rejuvenation

The GF2 Skin Rejuvenation formulation contains Hypochlorous acid aligned to the pH of the skin. In addition to providing protection against the Coronavirus, with ongoing use, the GF2 Skin Rejuvenation formulation works to reverse environmental damage such as UV exposure, pollution, toxins in many skincare products and normalise skin-cell functioning

### How to use GF2 Skin Rejuvenation

For protection against infection due to germs such as harmful bacteria and viruses, spray GF2 on the hands and face frequently when in public or travelling.

For best results:

- First clean the skin, removing any creams or makeup.
- Spray the skin with GF2. Let it dry naturally without rubbing or touching the treated area.

The active ingredient is completely natural and safe. It can be used:

- in the mouth
- on the skin
- in the eyes
- by pregnant and breastfeeding mothers
- on babies and children.

### Safety of HOCl when compared to other forms of chlorine

HOCl is classified as a Food Additive - USA FSIS (Food Safety Inspection Service). The HOCl molecule has been extensively researched and multiple peer-reviewed medical publications refer to its efficacy and safety, when used outside of the human body.

# INITIATIVE COVID-19

## Products Supply Emergency



	GF1	GF2	Trifectiv Plus
<b>Concentration HOCl</b>	350 mg/l	250 mg/l	370-400 mg/l
<b>pH</b>	5.5 – 6.0	5.5 – 6.0	5.5 – 6.0
<b>Volume</b>	50 ml, 500ml Professional	100 ml	100 ml, 500 ml Professional
<b>Class Medical Device</b>	C (FDA IIb) Although the product has medical device certification, we do not apply this in terms of the clinical use	A (FDA I) Although the product has medical device certification, we do not apply this in terms of the clinical use	C (FDA IIb) Certified medical device with SAHPRA, CE Mark and FDA certification planned for the next 12 months
<b>Application</b>	As a skin spray. For use on skin during and after (24 hours) aesthetic intervention. Prep of the skin as an antiseptic. Control of inflammation and to accelerate healing and reduce downtime. Prevention of PIH	As a skin spray. For use on intact skin as a maintenance treatment. Skin rejuvenation. Control of inflammation. Prevention of pigmentation disorders. Improvement of skin hydration and appearance.	Open wounds, wound irrigation, irrigation after elective surgery Lung nebulizing
<b>Provides extra HOCl for chlorine demand</b>	Yes	Yes	Yes
<b>Indications</b>	Disinfection of the skin before aesthetic intervention. Aftercare for skin manipulation e.g. laser, peelings, ablation treatment, needling, threads, filler injections, PRP	No organic bio matter present. For use in skin rejuvenation, radiation dermatitis, rosacea, seborrheic dermatitis (without open sores), acne, superficial skin infection - bacterial, viral and fungal. Sunburn. For the control of chronic low-grade inflammation.	Open wounds, wounds with bio load (plasma, puss, blood, organic matter). For nebulizing with a mesh nebulizer for lung infections. Nappy rash, umbilical cord infection, cradle cap. Severe impetigo. Seborrheic dermatitis (with open sores). Tonsillitis, aphthous ulcers mouth
<b>Contra indication</b>	Do not use with other "aftercare" products. Use on a clean skin. Do not use occlusive dressings or ointments with GF1	Do not mix with other products. Only apply other products after the skin has dried after application.	Uncontrolled diabetes, continuous presence of necrotic material in the wound, undebrided abscess. Do not use on the face
<b>Composition</b>	<i>50 ml</i> <i>500 ml</i>	<i>100 ml</i>	<i>100 ml</i> <i>500 ml</i>





### **CATALOG PRODUCTS COVID -19**

Work with Velas Largas means being able to access the consolidated international network consisting of industrial partnerships with high monthly production capacities, high quality products with certifications (CE and FDA) and valid Conformity Test Reports that offer peace of mind during the purchase phase not only on the product but also on the delicate phase of imports and acceptance of goods by different international customs.

Velas Largas also offers, through this partnership, a service

### **CONTACTS**

<https://www.velaslargas.com>

Rua Dr. António Martins, Número 44, 2º Esq,  
Benfica, Lisboa, Portugal

Tel.+39 3355767181  
info@velaslargas.com