



RayFocus

*Recover faster.
Perform better.*



Low Level Laser Massager


1. Modular: augmented use to all body parts
2. Beneficial: pain relief, anti-inflammation, metabolism enhancement
3. Convenient: wearable, rechargeable
4. Portable: compact, lightweight
5. Durable: used for more than 20 sessions after charging
6. Safe: patented laser safety design
7. Effective: red LED, red/NIR LD, massaged from shallow to deep
8. High-performance: 1199 Hz B2 Bahr frequency, PWM
9. Stylish: magnetic suction, easy to be stored in a personal bag



《Color codes》

Specifications

- Model: RayFocus1
- Light sources: 808nm LD 240mW, 660nm LD 100mW, Red LED 130mW
- Frequency: B2 Bahr frequency 1199Hz
- Laser safety: 3R
- Laser safety protection: YES (Laser emits only if skin touched)
- Power: Li polymer battery 1000mAh
- Charging port: USB-C
- Duration per session: 10min
- Sessions able to be used after fully charged: >20
- Weight: 42 (g)
- Size: 58*26*26 (mm)
- Content: Light engine *2, charging box *1, belt *1



Disclaimer

- Keep out of reach of children to prevent injuries to their eyes.
- Avoid looking directly at the light.
- Do not use it if you have a pacemaker, as the laser may cause interference with the performance of the pacemaker.
- Not intended to diagnose, treat, cure, or prevent any disease or health condition.